**REF** 193P, 194P, 195P

# Gelli-Roll®

## Operation Instructions

Models 193P, 194P, and 195P

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Model 193P, 194P en 195P

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Modell 193P, 194P und 195P

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Modelos 193P, 194P y 195P

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Modelli 193P, 194P e 195P

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Modelos 193P, 194P e 195P

## Betjeningsanvisninger

Model 193P, 194P og 195P

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Gentherm Medical, LLC  
12011 Mosteller Road  
Cincinnati, Ohio 45241, U.S.A.  
Tel: 1-800-989-7373  
[www.gentherm.com](http://www.gentherm.com)

Authorized European Representative  
CEpartner4U, BV  
Esdoornlaan 13  
3951 DB Maarn  
The Netherlands  
[www.CEpartner4U.com](http://www.CEpartner4U.com)

56289-N  
(06/2022)

CE 0344

## WARNINGS

- Physicians order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a physician. Higher risk patients; including pediatric, temperature sensitive patients and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to all bony prominences. Notify physician for changes in patient condition.
- Follow physician order/hospital policy for discontinuation of therapy. Notify physician for changes in patient condition.
- The area between the patient and the pad should be kept dry to avoid injury to the patient. Prep solutions have been reported to injure the skin when allowed to remain between patients and a water-circulating heating pad during prolonged procedures.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Do not use the device distal to arterial cross clamping. Non-observance can lead to thermal injury.
- Leaky pads or hoses should never be used. Water leaks present a risk of patient infection and should be handled accordingly.
- Do not place additional heat sources between the patient and the Gelli-Roll®. Tissue damage may result.
- Do not block fluid pathways with positioning accessories. Patient therapy may be interrupted.
- Do not position Gelli-Roll® at an incline. Potential for unexpected patient movement.
- Proper cleaning procedures must be maintained to prevent cross contamination.
- Do not place sharp objects in contact with pad as punctures may occur. Do not use pins or any instruments with sharp tips to secure pad or objects to the pad. Water leaks present a risk of patient infection and should be handled accordingly.

## PRECAUTIONS

- Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare professional.
- Presence of water or gel on the pad while the pad is not operating may cool the patient.
- Use of unapproved equipment or hoses is not recommended. The Gelli-Roll® has been designed and tested for use with GENTHERM thermal regulating systems only.
- Do not use alcohol to clean the Gelli-Roll®. Alcohol may cause pad and unit deterioration.
- Conductive pad therapy effectiveness will be reduced if pad is not in direct contact with the patient. Do not place multiple layers of sheets, or other objects that will reduce heat transfer between the pad and the patient.
- The Gelli-Roll® reusable hyper-hyperthermia blanket has a useful life of two (2) years.

## INDICATIONS FOR USE

- The Gelli-Roll® reusable hyper-hypothermia blanket is intended to be used by trained healthcare professionals as part of GENTHERM's thermal regulating systems, specifically the Blanketrol® II, Blanketrol® III, Hemotherm®, Norm-O-Temp®, or Micro-Temp® LT\* to lower or to raise a patient's temperature and/or maintain a desired patient temperature through conductive heat transfer.
- \*Note: Only the Cat#193P has been approved for use with the Micro-Temp® LT device.

## CONTRAINDICATIONS

- None Known.

## SET UP

1. Place the rolled Gelli-Roll® on the bed, stretcher or OR table. Unrolling the pad prior to the warming period may increase the time needed for the pad to reach the desired temperature.
2. Connect the male coupling and the female coupling from the Cat#286 hose to the unit. Make sure they “click” as they lock together.
3. Attach the pad to the Cat#286 hose by inserting the pad connectors into the hose connectors. Connectors should “click” as they lock together.
4. Check the equipment reservoir to ensure it is full of sterile water or water that has been passed through a filter less than or equal to 0.22 microns. Turn the unit on and select the desired temperature setting. The unit should be turned on approximately 45 minutes prior to first use, preferably first thing in the morning.
5. Setting of temperature is dependent upon patient’s core temperature and patient’s clinical condition. This temperature should be adjusted as required to attain and/or maintain target core temperature according to the physician’s order.
6. After the allotted time, unroll the Gelli-Roll®. Ensure that the pad lays flat with the hose routed without kinks towards the unit. Check pad before each use for product integrity (i.e. no leaks, flow is normal, etc.).
7. If cooling mode is selected, cool-down will take approximately 40 minutes.
8. After filling and positioning, make sure the pad and hose are free of kinks that may restrict flow.

## OPERATING INSTRUCTIONS

- Follow the thermal regulating system operation manual for operating instructions.
- To discontinue therapy: turn the regulating system off. The water in the pad will automatically drain back into the unit. Allow it to fully drain before disconnecting the pad from the unit or removing from bed.

## CLEANING

- For cleaning and disinfecting, always use conventional hospital-approved topical equipment cleaners and disinfectants that do not contain alcohol. Pads are washable at temperatures of approximately 108°F (42.2°C). Make sure pads do not enter washing machines or dryers. Avoid alcohol and other strong, undiluted disinfectants. These may cause staining or hardening of the pad’s outer skin. Thoroughly rinse product with clear water to remove any residue from cleaning solutions. A solution composed of 1:10 parts bleach to water can also be used for cleaning purposes.

**Note:** Do not use gas sterilization or autoclaving for cleaning and disinfecting this product.

## MAINTENANCE AND SERVICE

### Level I Repair

(Use GENTHERM Resurfacing Roll #169 included with pad)

Maintenance and service activities will sometimes overlap. In general, maintenance refers to any activity that does not require a certified technician. Maintenance may be performed by healthcare personnel or by other trained persons.

The following actions are considered maintenance:

1. Inspecting, cleaning, and disinfecting pads, and hoses.
2. Replacing hoses, other accessories.

Service refers to any activity that requires a Medical Equipment Service Technician, Certified Biomedical Electronics Technician, or a Certified Clinical Engineer.

The following actions are considered service:

1. Repair of Gel Pad.
2. Replacing hoses, and other accessories.
3. Other Repairs.

For repair of the Gel Pad Only. Using Resurfacing Roll provided:

Do not attempt to resurface any product that has multiple damaged areas, deep punctures, gouges, large cracks or tears, or areas where polymer has become dislodged. Items like this should be discarded and replaced.

1. Area to be resurfaced must be clean and dry before applying clear tape.
2. Any sticky surface exposed at damaged area may be removed with cleaning agent or warm, soapy water. Dry surface completely before applying clear tape.
3. Measure damaged area and cut tape to size. Allow 0.25" overlap on all sides.
4. Carefully peel back approximately one-half of brown release paper and crease it.
5. Carefully apply exposed portion of clear tape to damaged surface. Push out any air bubbles, that may be trapped beneath the tape.
6. Continue to carefully apply clear tape while slowly removing paper backing. Attempt to push out any air bubbles that develop.

### Other Repair

If a fold / kink is observed in the internal pad:

1. Lay the Gelli-Roll® flat on a solid surface.
2. Hold the Gelli-Roll® by the foot of the pad.
3. Gently shake the Gelli-Roll® like a small carpet until the fold / kink is removed.

## STORAGE

- Always handle Gelli-Roll® pads carefully. Never tear them away from the bed, stretcher, or OR table. Lifting the product from the corners may stretch the outer skin, causing stretch marks or tears. Pads should be rolled and cradled (like a baby) when lifted or carried. Carrying them unsupported can stretch the outer skin and cause undue stress that may result in splitting. Carts can be used for larger pads. These products are constructed to existing standards of quality control and should be usable for long periods with proper care.
- Pads should be stored in a dry area, in flat or rolled position (never folded). Keep pads out of direct sunlight. Keep cabinet doors and other sharp objects from coming in contact with pad surfaces at all times.

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