EL

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
Thermal Regulating System

WARNINGS

- Licensed healthcare practitioner's order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a physician. Higher risk patients; including pediatric, temperature sensitive and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to bony prominences. Notify licensed healthcare practitioner for changes in patient condition.
- Do not place patient on insulated tubing. Patient injury may occur.
- Follow licensed healthcare practitioner order/hospital policy for discontinuation of therapy. Notify licensed healthcare practitioner for changes in patient condition.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Leaky vests or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not block the fluid pathways with positioning accessories. Patient therapy may be interrupted.
- Do not place sharp objects in contact with the Maxi-Therm® Lite Patient Vest as punctures may occur. Do not use pins to secure the vest.
- The Maxi-Therm® Lite Patient Vest is single patient use only. Vests should be changed if it is used more than 7 days on the same patient or if it gets contaminated. Reuse may result in biocontamination or cross-infection.

PRECAUTIONS

- Caution: Federal law restricts the device to sale by or on the order of a licensed healthcare practitioner.
- Presence of water or gel on the vest while vest is not operating may cool the patient.
- Use of unapproved equipment or hoses is not recommended. The Maxi-Therm® Lite has been designed and tested for use with Gentherm thermal regulating systems only.
- Conductive pad therapy effectiveness will be reduced if the pad is not placed in direct contact with the patient.
- If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury.
- Do not place bean bags or other firm positioning devices under the blankets/pads. Skin damage may occur.
- In surgery, heating may affect toxicity of certain prep solutions which may cause skin injury when remaining on patients' skin during application of prolonged heat.
- Keep the area between the patient and the blanket dry. Excessive moisture may result in skin damage.

INDICATIONS FOR USE


The Maxi-Therm® Lite hyper-hypothermia vest is intended to be used by trained healthcare professionals as part of Gentherm's thermal regulating systems, specifically the Blanketrol® II or Blanketrol® III to lower or raise a patient's temperature and/or maintain a desired patient temperature through conductive heat transfer.

CONTRAINDICATIONS

None known.

SET UP

- Unfold vest on bed, stretcher, or operating room (OR) table with hose routed without kinks towards the unit. Check vest before each use for product

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