EL

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WARNINGS

- Licensed healthcare practitioner’s order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a licensed healthcare practitioner. Higher risk patients; including pediatric, temperature sensitive and operating room patients should be checked more frequently.
- Monitor patient skin temperature and condition paying particular attention to bony prominences. Notify licensed healthcare practitioner for changes in patient condition.
- Follow licensed healthcare practitioner order/hospital policy for discontinuation of therapy. Notify licensed healthcare practitioner for changes in patient condition.
- Do not use with a patient with an ischemic limb. Thermal injury may occur.
- Leaky blankets or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
- Do not place additional heat sources between the patient and PlastiPad®. Tissue damage may result.
- Do not block the fluid pathways with positioning accessories. Patient therapy may be interrupted.
- Do not place several layers of sheets or other materials that will reduce the heat transfer between the blanket and the patient.
- Do not place sharp objects in contact with the PlastiPad® blanket as punctures may occur. Do not use pins to secure PlastiPad® blanket.
- Do not use the device distal to arterial cross clamping. Non-observance can lead to thermal injury.

PRECAUTIONS

- Caution: Federal law restricts the device to sale by or on the order of a licensed healthcare practitioner.
- Presence of water or gel on the blanket while the blanket is not operating may cool the patient.
- Use of unapproved equipment or hoses is not recommended. The PlastiPad® has been designed and tested for use with Genthern thermal regulating systems only.
- If warranted, use pressure relieving or pressure reducing devices under the patient to protect from skin injury.
- Do not place bean bags or other firm positioning devices under the blankets/pads. Skin damage may occur.
- In surgery, heating may affect toxicity of certain prep solutions which may cause skin injury when remaining on patients’ skin during application of prolonged heat.
- Keep the area between the patient and the blanket dry. Excessive moisture may result in skin damage.
- The PlastiPad® reusable hyper-hypothermia blanket has a useful life of one (1) year.

CLEANING

Maintenance may be performed by nursing personnel or by other trained persons. The following actions are considered maintenance: Replacing hoses, cleaning hoses, and blankets.

- Clean blanket with a mild soap and water solution. Disinfect according to hospital/institution protocol. Genthern does not recommend sterilization. Do not autoclave.

STORAGE

- To store blanket, gently fold the corner containing the tubing towards the center of the blanket. Roll or fold blanket loosely avoiding sharp creases.

- Store blankets in cool, dry environment free from items that may puncture the blanket.

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