Aikuisten aluspeite

VAROITUKSIA:

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



FI

JA

					
Adult Underbody Blanket	Thermal Regulating System				
					

Adult Underbody Blanket

WARNINGS

- Physician’s order is required for temperature setting and for use of therapy. Patients should be checked at least every 20 minutes or as directed by a physician. Higher risk patients; including pediatric, temperature sensitive patients, and operating room patients should be checked more frequently.

- Failure to monitor the patient’s condition during therapy may result in skin damage or an undesirable elevated patient temperature.

- Follow physician order/hospital policy for discontinuation of therapy. Notify physician for changes in patient condition.

- Reduce air temperature or discontinue warming therapy per physician order. Notify the physician for changes in patient condition.

- Treatment of patient with poor circulation should be started in the lower temperature ranges (LOW or MEDIUM setting). The temperature setting can be increased as tolerated. In case of impaired circulation, use reduced temperature settings.

- Do not use the FilteredFlo® blankets distal to arterial cross clamping. Non-observance can lead to thermal injury.

- Do not use with a patient with an ischemic limb. Thermal injury may occur.

- FilteredFlo® Blankets are intended for single patient use only. Reuse may result in biocontamination or cross-infection.

- Use of Warming unit not connected to blanket may result in thermal injury.

PRECAUTIONS

- Caution: US Federal law restricts this device to sale by or on the order of a licensed healthcare professional.

INDICATIONS FOR USE

The FilteredFlo® Adult Underbody blanket is intended to be used with the WarmAir® patient warming system. The WarmAir® patient warming system is intended to prevent hypothermia and/or reduce cold discomfort before, during, and after surgical procedures. The thermal regulating system is used to raise a patient’s temperature and/or maintain a desired patient temperature through convective heat transfer from the controller to a warm-air heated blanket. The single-patient use blankets transfer the thermal energy to adult, pediatric, or neonate patients to obtain/maintain normal body temperature. It is intended for use by appropriately trained healthcare professionals in clinical environments.

CONTRAINDICATIONS

None known

PRODUCT DESCRIPTION

The Underbody blanket is part of the FilteredFlo® product line. It is positioned under the patient and is used with WarmAir® units for convective heat transfer.

SET UP

- Remove the Underbody Blanket from the package and place on the table.
- Unfold the Underbody Blanket, with the white material touching the table.
- Place patient on the side of the blanket that is labeled “this side to patient”. Insert WarmAir® hose into air inlet port.
- Patient’s head can be placed at either end of the blanket.

OPERATING INSTRUCTIONS

- Follow the WarmAir® Operation Manual for operating instructions.

CLEANING

- FilteredFlo® blankets are single-patient use and should be discarded after each use and per hospital protocol.

STORAGE

- Store blanket in cool, dry environment, free from items that may puncture the blanket.

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